



K&S Associates, Inc.
 1926 Elm Tree Drive
 Nashville, Tennessee 37210-3718
 800-522-2325 Fax 615-871-0856



**NON-INVASIVE KVP METER/MULTI-FUNCTION METER
 CALIBRATION SERVICE
 ORDER FORM**

FACILITY: _____
 CONTACT PERSON: _____
 INSTRUMENT MFGR / MODEL: _____
 FAX: _____ PHONE: _____

Calibration units and ranges available:

Diagnostic Radiology Range - W target, Al filter X-ray tube, 50 - 150 kVp, 3.4 mm Al HVL at 80 kVp
Mammography Range - Mo target, Mo filter X-ray tube, 22 - 35 kVp, 32 μ m, Mo filter

ORDERING AND PRICES - Please fill out one form for each kVp meter to be calibrated.

Standard Calibration: 6 kVp calibration points across one range...

Timer accuracy	_____
mA/mAs meter accuracy	_____
Dose (relative dose) linearity	_____
External chamber calibration (first point)	_____
additional calibration point	

Calibration intervals:

Your re-calibration reminder will be set up on a "manufacturer's recommended" interval unless you specify otherwise. IF you wish a re-calibration date recorded on your report please indicate the interval you require in the space provided.

 Phone 800-522-2325 Fax 615-871-0856
 email: ksinfo@kslab.com
 web: <http://www.kslab.com>
 for more information about K & S Associates and our services
 please call or visit our web site